



# Nucleic Acid Extraction

# KingFisher product portfolio

## KingFisher Duo Prime



12 - 24 samples 96 well plates  
6-12 samples 24 well plates  
30-5000 µl

## KingFisher Flex



96 samples  
24 samples  
20-5000 µl

## KingFisher Presto



24/96 samples  
Full Automation with robots  
50-5000 µl

## Plastic consumables



## BindIt Software



## KingFisher & MagMAX kits



# Versatility Through Various Throughput Levels

KingFisher instruments provide throughput options that suit a wide variety of research needs

## Automation “lite”



## Semi-automation



## Full automation



### Duo Prime system

### Flex system

### Presto system

	Duo Prime system	Flex system	Presto system
<b>Instrument size</b>	Stand-alone compact benchtop	Stand-alone benchtop	Benchtop + robotic liquid handler
<b>Throughput level</b>	● Low to medium	● ● Medium to high	● ● ● Ultrahigh
<b>Samples per run/plate format</b>	12 or 6 per run	96 or 24 per run	96 or 24 per run



# Reagents and Kits Across a Broad Spectrum of Application Types

DNA	RNA	Total nucleic acid	Immunoprecipitation / Protein
gDNA	Total RNA (excluding small RNA)	Viral RNA/DNA	Protein A, G
Cell-free DNA	True total RNA (miRNA and large)	FFPE DNA/RNA	Secondary antibodies
	mRNA	Cell-free total nucleic acid	Streptavidin

Find the right kit at [thermofisher.com/kingfisherkits](https://www.thermofisher.com/kingfisherkits)

# KingFisher Consumables

- Specially designed KingFisher microplates, tubes and tipcombs for optimal processing and isolation of target molecules
- Made of polypropylene
- Optimal for magnetic particle processing due to low binding of biomolecules
- Enables >99% (KingFisher) and >95% (KF mL & KingFisher Flex) recovery of magnetic particles





## Applications for molecular biology

# KingFisher Applications

## ***DNA***

gDNA  
Viral DNA  
Plasmid DNA  
PCR Clean up  
Phage Display

## ***Protein***

Titanium beads  
Streptavidin-Biotin  
Thermo Pierce Beads  
**Trypsin SMART Digest**

## ***RNA***

Total RNA  
miRNA  
Poly A Selected



## ***Mag Beads***

Luminex  
Omega Biotek  
Spherotech Beads  
Sera Mag  
Mag sepharose  
Promega  
Perkin Elmer


## ***Dyna Beads***

Immuno Percipitation (IP)  
Cell Isolation (+ve and -ve)  
T-cell activation beads  
Custom Beads  
IVD Assay Development

## ***Instruments***

NGS - (IonTorrent and MiSeq)  
CTC - IsoFlux/Liquid Biopsy  
Chromatography - Dyna RPC 18  
BioRad – Bio-plex  
Mass Spec - Upstream isolation,  
clean up or digest


# Nucleic Acid Extraction Process






# Qubit assays

- The Qubit uses target-selective dyes that emit fluorescence when bound to DNA, RNA or protein
  - The system won't detect other biomolecules in the sample
- DNA quantification assays: dsDNA, ssDNA and oligos
- RNA and miRNA quantification assays
- Protein quantification assays



# Qubit assays

- All assays use the same general protocol
- 2min incubation time for DNA and RNA assays
- 15min incubation time for protein assays



## Qubit 4



## Qubit 4 Quantification Starter Kit




## Qubit 4 NGS Starter Kit




# ThermoFisher SCIENTIFIC

---





# Molecular Biology Brand Strategy: reagents



- **Portfolios:**

- End Point PCR (ePCR) – includes enzymes and Master Mixes for end point PCR and reverse transcription, and nucleotides
- Thermal cyclers and PCR plastics
- Molecular Biology Tools (MBT) – includes enzymes, products for cloning and NA electrophoresis, competent cells



# Combined Portfolio Expands Addressable Markets



# PCR is everywhere

96% of all mol bio labs are performing PCR



A step critical for most mol bio applications



PCR is a significant sales opportunity



## Milions

PCR reactions performed in Africa every year



Cloning/gene characterization



Genotyping



Sequence detection



Sequencing



# PCR components



**Instruments**



**PCR Plastics**




**PCR enzymes**



**dNTPs**

# Reverse transcription

## Invitrogen™ Superscript™ IV Reverse Transcriptase



- Super efficient - high cDNA yields and processivity
- Super robust - transcribes even from degraded or inhibitor-containing RNA samples
- Super convenient - available in a variety of formats, all with a fast, 10-minute cDNA synthesis protocol

[thermofisher.com/ssiv](https://thermofisher.com/ssiv)

## Invitrogen™ Superscript™ IV VILO™ Master Mix



- Improved Ct—up to 8 fewer cycles compared to other RT reagents
- Superior linearity—works with a broad range of input RNA for any RT-qPCR application

[thermofisher.com/4vilo](https://thermofisher.com/4vilo)

## SuperScript™ IV One-Step RT-PCR System



- Two-phase hot-start activation mechanism—enables room temperature setup and high specificity
- Unmatched performance—high sensitivity and target length, with the fastest protocol

# DNA Polymerases

## Platinum II Taq Hot-Start DNA Polymerase



- **Universal primer annealing—reduces tedious optimization**
- **Engineered Taq polymerase combined with leading hot start—enables fast, robust, and highly specific PCR**

[thermofisher.com/platinumiiatq](https://thermofisher.com/platinumiiatq)

## Invitrogen Platinum SuperFi DNA Polymerase



- **Exceptional fidelity—100x more accurate than Taq DNA polymerase**
- **Robust and versatile—ideal for difficult targets**

[thermofisher.com/superfi](https://thermofisher.com/superfi)

## Thermo Scientific DreamTaq Hot Start DNA Polymerase



- **Better reaction outcomes**
  - Higher yields
  - Higher specificity and sensitivity
  - Longer amplicons supported
- **Convenience**
  - Room temperature reaction setup
  - Direct loading on a gel with green buffer formats
  - Minimal reaction optimization due to optimized buffer with MgCl<sub>2</sub>



**DreamTaq Hot Start enhances your PCR productivity**

# Thermo Scientific Phusion High-Fidelity DNA Polymerase

## Highlights:

- **Accurate sequence**—high fidelity (52x more accurate than Taq DNA polymerase, 6x more accurate than Pfu DNA polymerase)
- **Reputation**—most widely used high-fidelity enzyme on the market
- **Formats optimized for specific applications**—Phusion U DNA Polymerase for dUTP incorporation, green buffer for direct gel loading, and multiplex master mix formats available
- **Money-saving**—affordable for a high-fidelity DNA polymerase

[thermofisher.com/phusion](https://thermofisher.com/phusion)



# Analyse

**Christian Dutoziet**  
TSS – West and Central Africa

1

Introducing E-Gel Power Snap Electrophoresis System

2

Comparison

3

Consumables



**Instruments**



**Precast Gels**



**Ladders and Stain**



**Plastics**

# Introducing E-Gel Power Snap Electrophoresis System



The E-Gel™ Power Snap Electrophoresis System provides:

- **Faster analysis**—sample loading to image capture in as little as 15 minutes
- **Safer Choice**—minimize handling of hazardous chemicals as well as risk of gel breaking
- **Simple operation**—intuitive user interface with large touch screen and integrated operating system
- **Streamlined electrophoresis system**—with integrated image visualization
- **Small footprint**—for convenient benchtop use

# E-Gel Power Snap Electrophoresis System

The only integrated gel running and imaging platform

**High resolution digital camera**  
for quick, easy image capture

**Touch screen**  
with intuitive user interface

**Intuitive analysis software**  
with auto-exposure and auto-capture

**Front facing USB port**  
for easy file transfer

**32 Gb internal storage**  
>3 months of work with camera

Designed for standard  
**E-Gel Cassette electrophoresis**

Compatible with all standard low-throughput E-Gel Agarose gels, including: E-Gel EX, E-Gel SYBR Safe, E-Gel SizeSelect and E-Gel CloneWell, E-Gel Go!




**Amber filter**  
for real-time sample monitoring

Integrated power unit and  
**blue light transilluminator**

**Pre-programmed protocols**  
for rapid analysis

**E-Gel Go! Adapter**  
to enable low throughput E-Gel Go! Agarose gels





# Questions?

